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aeruginosa carriage in residents of French and German long-term care facilities. *Clinical Microbiology and Infection*, 28(10), 1353–1358. <https://doi.org/10.1016/j.cmi.2022.05.004>

3. Tu, M. G., Lin, C. C., Chiang, Y. T., Zhou, Z. L., Hsieh, L. Y., Chen, K. T., Chen, Y. Z., Cheng, W. C., & Lo, H. J. (2022). Distribution of Yeast Species and Risk Factors of Oral Colonization after Oral-Care Education among the Residents of Nursing Homes. *Journal of Fungi*, 8(3). <https://doi.org/10.3390/jof8030310>
4. Weppner, J., Gabet, J., Linsenmeyer, M., Yassin, M., & Galang, G. (2021). Clostridium difficile Infection Reservoirs Within an Acute Rehabilitation Environment. *American Journal of Physical Medicine & Rehabilitation*, 100(1), 44–47. <https://doi.org/10.1097/PHM.0000000000001579>
5. Rowan, N. J., Kremer, T., & McDonnell, G. (2023). A review of Spaulding's classification system for effective cleaning, disinfection and sterilization of reusable medical devices: Viewed through a modern-day lens that will inform and enable future sustainability. *Science of The Total Environment*, 878, 1–18. <https://doi.org/10.1016/j.scitotenv.2023.162976>
6. Spaulding, E. (1970). The role of chemical disinfection in the prevention of nosocomial infections. *International Conference on Nosocomial Infections*.
7. Klacik, S. (2019). 'Noncritical' OR items deserve due diligence in disinfection - OR Manager. *International Association of Healthcare Central Service Materiel Management*. <https://www.ormanager.com/noncritical-items-deserve-due-diligence-disinfection/>
8. Ontario Agency for Health Protection and Promotion (Public Health Ontario), & Provincial Infectious Diseases Advisory Committee. (2013). *Best practices for cleaning, disinfection and sterilization of medical equipment/devices in all health care settings* (3rd ed.). Queen's Printer for Ontario.
9. Werkgroep Infectie Preventie, de Bruyn, A. C. P., van Klingeren, B., & Severin, W. P. J. (2009). *WIP richtlijn - Beleid reiniging desinfectie en sterilisatie*. www.wip.nl
10. Gardner, P., Muller, M. P., Prior, B., So, K., Tooze, J., Eum, L., & Kachur, O. (2014). Wheelchair cleaning and disinfection in Canadian health care facilities: 'That's wheelie gross!' *American Journal of Infection Control*, 42(11), 1173–1177. <https://doi.org/10.1016/j.ajic.2014.08.007>

11. *Cushioning Products by Supracor - User manual - Stimulite*. (2018). www.supracor.com
12. *Jay Foam Cushion Products User Instruction Manual & Warranty*. (2021).
13. *JAY J3 Cushion Owner's Manual*. (2017).
14. *Jay Xtreme Active Owner's Manual*. (2018).
15. *ROHO® DRY FLOATATION® Wheelchair Cushions Operation Manual*. (2019).
16. *Vicair Vector O2 User Manual*. (2018).
17. Op De Beeck, J. (2023). *Test Report Disinfection Wheelchair Cushion*.
18. Mohapatra, S. (2017). Sterilization and Disinfection. In *Essentials of Neuroanesthesia* (pp. 929–944). Elsevier. <https://doi.org/10.1016/B978-0-12-805299-0.00059-2>

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