alber®

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e-motion®

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e-motion®

Information for Therapists and DME Dealers

Table of Contents

1. General instructions	13
2. Adjustment options	14
3. Moving off	18
4. Individual adaptation of power assistance	20
5. Compensation for uneven application of force	20
6. Important instructions relating to range	21

Instructions on adaptation for specialised dealers and therapists

1. General instructions

Only optimum adaptation can lead to the achievement of the greatest possible freedom of movement for the user when it comes to actively power-assisted wheelchairs. e-motion provides a variety of options to help achieve harmonious adjustment in this respect.


The following functions and characteristics enable you to individually adapt to the wheelchair, whether employed individually or in combination:

Variably-mountable wheels


The fully-integrated arrangement of the motors and function elements in the wheel hub allow the realisation of variable wheel positions, depending on the wheelchair frame involved.

For example, a seat with a greater slope can be provided in the initial stages of mobilising the patient. On the other hand, it is also possible to react to different body positions. The latter is critically important if efficient propulsion and a good seating position is to be achieved.


Exampels:

Greater angle of seat
Mobilisation

Wheel position at front
Manoeuvrability


Low seat height
Movement using the tips of the toes

You can utilise the complete „active“ wheelchair adaptation spectrum and still have freedom of choice when it comes to selecting the most suitable wheelchair model.

The following sections are intended as a step-by-step introduction to optimum adaptation of e-motion.

Caution:

Please read the user instructions first and familiarise yourself with the function elements and safety instructions!

2. Adjustment options


Power assistance

The motor output can be adjusted in two stages at the touch of a button [1]. This enables adaptation of power assistance to suit the individual requirements of the wheelchair user.

Assistance level 1 is usual for use indoors and normal operation outdoors.

Assistance level 2 provides additional motor power for a variety of manoeuvres (e.g. overcoming gradients, kerbstones and ramps).

See chapter 7 of the user instructions for more information on this topic.

Sensor sensitivity

Sensitivity can be infinitely adjusted to meet requirements using an adjusting wheel [2], individual adjustment of the left or right wheel being possible. This option enables the user to compensate for discrepancies in the drive force.

Works adjustment allows for maximum deflection of the sensor [3] through the application of minimal force on the push rims. Rotating the adjusting wheel [2] in the „+“ direction increases the sensitivity of the sensor, and the user requires a lower application of force to achieve movement.

Rotating in the „-“ direction reduces sensor sensitivity, and the user must apply greater force to move the wheelchair.


Sensor sensitivity adjustment options can be deactivated by loosening the grub screw [4] in the adjusting wheel. After loosening the sensitivity of the sensor is no longer altered by turning the adjusting wheel.

Adjusting power assistance [1] generally has no effect on the sensor setting. Re-adjustment of the sensor after altering power assistance and/or deactivating the drive is therefore unnecessary.

The following diagram illustrates clearly the ratio between the power required and the sensor setting:

- Minimum = Sensor adjusting wheel (2) is entirely within the „+“ range
- Average = Sensor adjusting wheel is in the middle of the adjusting range
- - - - Maximum = Adjusting wheel is turned fully into the „-“ range

Sensor sensitivity


Minimal residual force is sufficient in stage 2 (up to 100% motor assistance). A person weighing approx. 80 kg need only apply a force of approx. 500 gram (5 Newton's) on a level surface to effect movement from a stationary position. Even less force is subsequently required to maintain a continuous forwards or backwards movement.

Conversely, a force of approx. 3 kg (30 Newton's) is necessary to realise a movement when the adjusting wheel has been tightened completely (entirely within the „-“ range).

Individual sensor sensitivity selection also assists therapeutic measures. The user can be gradually conditioned to exert a greater force. This is achieved as a result of reduction in sensor sensitivity accomplished at regular intervals over a given period of time.

Works adjustment of wheels

e-motion wheels are adjusted by the manufacturer, the unit being delivered with power assistance set at stage 1 and a high sensor sensitivity level [3]. The adjusting wheel [2] is almost completely within the „+“ range.

Programming the drive parameters

The programming switch is located behind a sliding foil [1] within the wheel hub.

Rotary switch [2]

10 different recorded drive parameters (combinations of starting, degree of assistance and after-running parameters) can be adjusted with this feature, optimally adapting e-motion to suit the wheelchair user's requirements.

Drive parameters

Switch setting	Drive parameter (stage 1)			Drive parameter (stage 2)			
	Start*	Assistance*	Afterrunning*	Start*	Assistance *	Afterrunning*	
Works setting	0	0,3 sec	50	1,0 sec	0,3 sec	80	3,3 sec
	1	0,4 sec	40	0,8 sec	0,2 sec	70	4,2 sec
	2	0,5 sec	50	1,2 sec	0,3 sec	80	3,6 sec
	3	0,6 sec	60	1,2 sec	0,3 sec	80	3,4 sec
	4	0,7 sec	70	1,0 sec	0,4 sec	90	3,2 sec
	5	0,8 sec	70	1,0 sec	0,4 sec	90	3,0 sec
e.g. for quadriplegics	6	0,9 sec	80	0,8 sec	0,8 sec	100	2,8 sec
	7	1,0 sec	80	0,8 sec	0,8 sec	100	2,6 sec
e.g. for older patients	8	1,1 sec	90	0,6 sec	1,0 sec	100	2,4 sec
	9	0,3 sec	50	0,6 sec	0,3 sec	100	3,3 sec

* Explanations

Start:

The time period in seconds between moving the push rim and the motor engaging. The longer the period, the slower the motor reacts to the impulse from the push rims. Higher adjustment values are ideal for patients who, although physically strong, are limited with regard to co-ordination of their arms (e.g. spastics).

Support:

What is the maximum percentage of the potential motor output which can be utilised?

Max. torque: 30 Nm (=100%).


Greater assistance = reduced range,


lower assistance = greater range.

After-running:

How long does the motor „run on“ after actuating the push rims and releasing them again?

Obstacles (e.g. kerbstones) can be more easily overcome, and less propulsion is necessary when driving uphill.

[3] DIP switch

This switch can be used to pre-adjust several system functions.

Switch 1

The drive stage/power assistance can no longer be adjusted if the DIP switch is set at OFF. This can be a practical measure if the range of functions needs to be reduced in as far as possible for safety reasons (e.g. when dealing with older patients).

ON: Switch activated (works setting)

OFF: Switch deactivated

Switch 2

The required drive stage (1 or 2) can be specified with this switch if DIP switch 1 (see above) is set at OFF.

ON: Assistance level 2 activated (works setting)

OFF: Assistance level 1 activated

Switch 3

Any signals which are not of critical importance when it comes to safety can be deactivated at the user's discretion by setting the DIP switch at OFF. Error signals and warnings relating to flat batteries continue to be indicated by a beeping tone.

ON: Signals on (works setting)

OFF: Signals off

Switch 4

The deactivation time can be increased from 30 minutes to 120 minutes. The longer the deactivation time, the more current is consumed during shutdown and the more limited the overall range.

ON: 30 min (works setting)

OFF: 120 min

Switch 5 + 6 do not have any function at present (redundant/reserve)

Accessories

e-motion can be ideally adapted to suit the user's requirements through the use of practical accessories, such as spoke covering (a protective feature), rubber push rims (improved grip) or overturning supports (as a safety measure).

3. Moving off

Ensure that overturning supports are fitted to the wheelchair before moving off. This is necessary for reasons of safety. The user should then be introduced to handling e-motion with works settings (stage 1, „high“ sensitivity). Select level ground with a smooth surface to begin (e.g. stone, PVC or parquet flooring), avoiding carpet¹ if possible. It may be necessary to begin at stage 2 in the case of users with more limited strength.

¹ Carpet fibres predetermine a particular direction, usually leading to the wheelchair driving in a slight curve. The curve is exaggerated in the case of deep pile carpets.

Allow the user to drive a few metres forward, stop and then accelerate again.

Important!

Ensure that equal pressure is applied to the left and right push rims when traveling forwards in a straight line. The wheelchair user may find this difficult initially, particularly if a transition has been made from a manual wheelchair to an actively power-assisted chair.


The following points should be observed to exploit the product's characteristics fully:

- The push rim need only be actuated relatively briefly and lightly, as the drive mechanism does any work needed
- Do not „tug“ or guide the push rim too long
- Actuate the left and right push rims evenly and simultaneously when driving straight ahead

The following exercises help in the initial stages when learning the new method of driving:


1st exercise:

Driving straight ahead: push the push rims simultaneously with both hands. Push and stop.


Driving in a curve and turning on the spot should be practised when braking and delaying have been mastered.

2nd exercise:
Curves/Turning left and right

3rd exercise:
Turning on the spot

The system gives electro-motoric assistance to the braking procedure when manually delaying the push rim. This means that the occupant need not apply as much force when braking the wheelchair. Braking support also increases as the power assistance setting is raised (stage 2).

It is recommended that power assistance be set at stage „2“ if a longer incline is to be tackled. The user only requires a fraction of the force to stop the wheelchair, and the palms of the hands are not subject to as much heat during the manoeuvre.


4. Individual adaptation of power assistance

Power assistance can be increased from stage 1 to stage 2 after simple driving manoeuvres have been mastered. The system now reacts with considerably increased assistance than was the case in stage 1.

Increasing power assistance leads to a parallel increase in programmed „after-running“. After-running is the time period in which the motor continues to provide both electrical and motoric assistance and drive after a propulsion movement has been realised. The longer the motor's after running phase after each propulsion movement, the less often the user needs to push to cover a given distance. This effect is of particular advantage when dealing with gradients. The user need not push as often (in contrast to manual wheelchairs). This encourages a more fluid movement and saves energy reserves.

Evenly-balanced pushing is extremely important in stage 2 if an exactly straight forward movement is to be achieved. Wheelchair users often drive in a meandering fashion initially in stage 2. This is a result of both wheel after-running, as described above, and the individual movement pattern of the user. User tests have, however, shown that beginners can generally safely master the functions after 1 - 2 days.

Now let the user carry out the exercises described above with power assistance in stage 2.

5. Compensation for uneven application of force

The adjusting wheels ([2], see page 4) on both e-motion wheels should be synchronised to achieve an even driving movement, any differences in the force applied by the user being thus compensated for. This can be tested by driving over a straight stretch of ground. Movement of the wheelchair which is uneven, requiring more pushing on one side, indicates that a greater exertion of force is necessary on this side.

Increase or decrease the sensitivity of one or the other of the wheels to compensate for a discrepancy. The user should cover a short, straight distance after every adjustment to evaluate the effect achieved.

The „R“ and „L“ stickers supplied with the chair should be used if a different sensor setting is selected, applying them to the left and right wheels respectively. This helps avoid incorrect switching of the wheels when removing and remounted. e-motion wheels can be fitted to either side if the new settings selected are the same.

6. Important instructions relating to range

The greatest range can be achieved when driving at assistance level 1. Test drives under these conditions resulted in the achievement of a range of more than 12 km. Frequent acceleration and delays, along with continual steering correction, can considerably reduce the range. Bumpy streets, wind conditions, hilly surfaces and, most particularly, manoeuvring at assistance level 2 in confined areas also result in increased power consumption and a reduction in range.

Power consumption and, consequently, the range also depend considerably on the driving style of the user. Some wheelchair users propel themselves through a lot of short, brief pushing movements, during which they are continually running the push rings through their hands, an action which actually applies a braking action to the rims. The resulting sensor deflection against the direction of movement directs braking current into the motor. The wheelchair user is thus continually pushing and braking without necessarily noticing this. However, the range is considerably reduced by this style of driving.

This can now and again lead to ranges of less than 3 km being achieved when driving in the highest assistance level at maximum sensitivity. Sensor sensitivity for this user group should be reduced considerably in the „-“ direction using the adjusting wheel [2] (see page 4).

Roll resistance also has a considerable effect (e.g. when the air pressure in the wheelchair steering wheel tyres is inadequate). The range is halved if, for example, roll resistance is 6% instead of 3%. The range is also halved if, in addition to a roll resistance of 3%, a gradient of 3% is to be covered. A gradient of 12% means that the range is only 1/5 of that achievable over level ground.

It goes without saying that driving under practical conditions will always involve a mixture of gradients, inclines and level stretches of ground, the combination of which in turn influences the range achieved.

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Ulrich Alber GmbH