

On the basis of the tests undertaken, the sample(s) of the below product have been found to comply with the requirements of the referenced specifications at the time the test were carried out.

**Applicant Name & Address** : Invacare Corporation  
1 Invacare Way  
Elyria, OH 44035 USA

**Product Description** : Non-rigid body support units

**Relevant Standards** : ISO 10535:2006

**Acceptance Level** : No part shall show any sign of damage or wear.  
Any dimensional change shall be determined by measurement

**Results** : Pass

**Verification Issuing Office Name & Address** : Intertek-Grand Rapids  
4700 Broadmoor Ave. SE Suite 200  
Kentwood, MI 49512 USA

**Date of Test** : September 11, 2013 to June 6, 2014

**Report Number(s)** : 101301824GRR-001A  
101301824GRR-001B  
101301824GRR-001C  
101301824GRR-001D

**Products Tested:**

Non-rigid body support unit type:	Universal Low,
	Universal Standard,
	Universal High,
	Universal High Plus
Size:	XS, S, M, L, XL
SWL:	200 kg, 250 kg for XL
Maximum Washing Temperature:	95°C
Material:	100% Polyester

Non-rigid body support unit type:	Easy-Fit
Size:	XS, S, M, L, XL
SWL:	200 kg, 250 kg for XL
Maximum Washing Temperature:	95°C
Material:	100% Polyester

Non-rigid body support unit type	Stand Assist
	Standing Transfer Vest
	Transfer Stand Assist
Size:	S, M, L, XL
SWL:	200 kg
Maximum Washing Temperature:	95°C
Material:	100% Polyester

Non-rigid body support unit type:	Dress Toileting
	Dress Toileting High
Size:	S, M, L, XL
SWL:	200 kg, 250 kg for XL
Maximum Washing Temperature:	95°C
Material:	100% Polyester

Non-rigid body support unit type:	Comfort,
	Comfort Toileting
	Comfort Standard
Size:	XS, S, M, L, XL
SWL:	200 kg, 250 kg for XL
Maximum Washing Temperature:	95°C
Material:	100% Polyester

Non-rigid body support unit type:	Amputee
Size:	S, M, L, XL
SWL:	200 kg, 250 kg for XL
Maximum Washing Temperature:	95°C
Material:	100% Polyester

Non-rigid body support unit type:	Lifting Band
Size:	Short or Long
SWL:	200 kg
Maximum Washing Temperature:	95°C
Material:	100% Polyester

Non-rigid body support unit type:	Stretcher Sling
Size:	One Size Only
SWL:	200 kg
Maximum Washing Temperature:	95°C
Material:	100% Polyester

Intertek Testing Services NA, In. Grand Rapids (Kentwood) facility is accredited in accordance with ISO/IEC 17025 by the American Association of Laboratory Accreditation (A2LA) [Certificate No. 0078.02](#)

This Verification is for the exclusive use of Intertek's Client and is provided pursuant to the agreement between Intertek and its Client. Intertek's responsibility and liability are limited to the terms and conditions of the agreement. Intertek assumes no liability to any party, other than to the Client in accordance with the agreement, for any loss, expense or damage occasioned by the use of this Verification. Only the Client is authorized to copy or distribute this Verification. Any use of the Intertek name or one of its marks for the sale or advertisement of the tested material, product or service must first be approved in writing by Intertek. The observations and test results referenced from this Verification are relevant only to the sample tested. This Verification by itself does not imply that the material, product, or service is or has ever been under an Intertek certification program.

**NOTE : This verification is part of the full test report(s) and should be read in conjunction with it.**



---

**Signature**

**Name: Anthony Serge**

**Position: Department Manager**

**Date: June 9, 2014.**